

A Gesture-based Tool for Sterile Browsing of Radiology Images

Team ID - PNT2022TMID25935

Run The Application

At last, we will run our flask application

```
if __name__ == "__main__":  
    # running the app  
    app.run(debug=False)
```

Run The app in local browser

- Open anaconda prompt from the start menu
- Navigate to the folder where your python script is.
- Now type “python app.py” command

Navigate to the localhost where you can view your web page

```
(base) E:\>cd E:\PROJECTS\number-sign-recognition\Flask  
(base) E:\PROJECTS\number-sign-recognition\Flask>python app.py
```

Then it will run on localhost:5000

```
* Serving Flask app "app" (lazy loading)  
* Environment: production  
  WARNING: This is a development server. Do not use it in a production deployment.  
  Use a production WSGI server instead.  
* Debug mode: off  
* Running on http://127.0.0.1:5000/ (Press CTRL+C to quit)
```

Navigate to the localhost (<http://127.0.0.1:5000/>) where you can view your web page.